

2004

# Detection of pesticide residues in aquaculture products

Borlongan, Ilda G.

Aquaculture Department, Southeast Asian Fisheries Development Center

---

Borlongan, I. G. (2004). Detection of pesticide residues in aquaculture products. In Laboratory manual of standardized methods for the analysis of pesticide and antibiotic residue in aquaculture products. (pp. 23–46). Tigbauan, Iloilo, Philippines: Aquaculture Department, Southeast Asian Fisheries Development Center.

---

<http://hdl.handle.net/10862/1698>

---

*Downloaded from <http://repository.seafdec.org.ph>, SEAFDEC/AQD's Institutional Repository*

---

---

## **CHAPTER 2**

---

---

# **Detection of Pesticide Residues in Aquaculture Products**

*Prepared by*

Ilda G. Borlongan, Ph. D.

Scientist II

Aquaculture Department

Southeast Asian Fisheries Development Center

Philippines

## INTRODUCTION

Agricultural chemicals such as fertilizers and pesticides have made an important contribution to agriculture. Pesticides protect crops from pests and diseases. They have brought about large yield increases, and have helped ensure that the rise in food production has kept well ahead of the rise in population. However, there is a growing concern about the safe use of these chemicals, and the potential dangers to farmers who use them, the environment, and consumers. There is particular concern about pesticides, since almost all chemicals that can kill pests are also potentially damaging to human health.


Legislation requires that pesticide use is appropriately controlled and maximum residue levels (MRLs) not be exceeded. The level of pesticide residues in food raw materials is a measurable standard. But while residue analysis is essential for companies wishing to assure themselves that their products have been produced in accordance with best practice and within the law, it can be used to greatest effect when targeted at samples most likely to contain residues.

Reliable residue analytical methods are necessary to measure the magnitude of residue in a seafood, and to enforce legal residue limits (tolerances). Sample preparation and extraction, clean up of extracts and pesticide detection are the main procedures in pesticide residue analysis. There is an interplay among these factors which should be considered in the choice of a particular method.

# Determination of Pesticides Residue in Aquaculture Products

## Preparation of Samples

The fish sample is scaled (if scaly) and the head removed. Using a stainless knife, the flesh and other edible portions are removed from the bone and entrails. These edible portions are cut into small cubes or pieces and frozen until analyses.



Head and shells of tiger shrimp samples are removed and edible portions separated. These edible portions are cut into small cubes or pieces and frozen until analysed.



Mollusks (oysters, clams and mussels) are shucked or de-shelled (snails) and meat are separated and cut into small pieces and frozen until analysed.



Oysters



Snails



Clams

Seaweed samples are chopped using a knife or cut into small pieces using a pair of scissors and frozen until analysed.



The frozen samples are allowed to thaw under room temperature and then ground in a meat chopper. The ground samples are thoroughly mixed and representative samples are taken for analysis.

## Multi-residue Method

### Principle

The sample is blended with anhydrous sodium sulfate to disintegrate the solid and for it to combine with water in the sample. Fat is extracted from sample using petroleum ether. Extracts are purified by petroleum ether-acetonitrile partitioning, followed by chromatography on Florisil column. For some fractions, further clean up in Magnesia column may sometimes be needed. Pesticide residues in concentrated eluates are measured by gas chromatography.

### Reagents

- a.) Acetonitrile ( $\text{CH}_3\text{CN}$ ) - Some lots of reagent grade  $\text{CH}_3\text{CN}$  are impure and require distillation. Generally vapors from such lots will turn moistened red litmus paper blue when held over the mouth of a storage container. Purify  $\text{CH}_3\text{CN}$  as follows: To 4 L  $\text{CH}_3\text{CN}$ , add 1 mL phosphoric acid ( $\text{H}_3\text{PO}_4$ ), 30 g phosphorus pentoxide ( $\text{P}_2\text{O}_5$ ) and boiling chips, and distill in all glass apparatus at  $81-82^\circ\text{C}$ . Do not exceed  $82^\circ\text{C}$ .
- b.) Acetonitrile saturated with petroleum ether - Saturate  $\text{CH}_3\text{CN}$  with redistilled petroleum ether
- c.) Methanol - UCP, reagent grade or ACS
- d.) Alcoholic alkali solution (2%) - Dissolve 2 g potassium hydroxide (KOH) in methanol, and dilute to 100 mL.
- e.) Eluting solvent (6%)- Dilute 60 mL ethyl ether to 1 L with redistilled petroleum ether.
- f.) Eluting solvent (15%)- Dilute 150 mL ethyl ether to 1 L with redistilled petroleum ether.
- g.) Eluting solvent (50%)- Dilute 500 mL ethyl ether to 1 L with redistilled petroleum ether.


- h.) Ethyl ether- Redistilled at 34-35° C, and stored under nitrogen (N<sub>2</sub>).  
Must be peroxide-free by test.
- i.) Florisil, 60/100 PR grade, activated at 675° C, available from Floridin Co., 3 Pennsylvania Center, Pittsburgh, PA 15235. When 675 °C activated Florisil is obtained in bulk, transfer immediately after opening to 500 mL glass jars or bottles, with glass-stoppered or foil-lined, screw-top lids, and store in the dark. Heat for 5 h or more at 130° C before use. Store at 130° C in glass stoppered bottles or in air-tight dessicator at room temperature and reheat at 130°C after 2 days.
- j.) Hexane-Reagent grade, redistilled in all glass distillation apparatus
- k.) Magnesium oxide (MgO) - Adsorptive magnesia (Fisher Scientific Co. No. S-120). Treat as follows: Mix 500 g MgO with water to form a slurry. Heat on steam bath for about 30 min, and filter with suction. Dry overnight at 105-130° C and pulverize to pass No. 60 sieve. Store in jar with cover.
- l.) Magnesia-Celite mixture- Mix treated MgO with Celite 545 (1:1, w:w).
- m.) Petroleum ether- Reagent grade, redistilled in all glass distillation apparatus at 30-60° C.
- n.) Sodium sulfate- Anhydrous, granular

## Extraction

1. Weigh 25-50 g of thoroughly ground and mixed sample in a homogenizing beaker.







2. Homogenize the sample in high speed blender or homogenizer until thoroughly mixed.

3. Add 100 g anhydrous sodium sulfate ( $\text{Na}_2\text{SO}_4$ ) to combine with water and to disintegrate sample



4. Mix with a spatula and blender until well mixed. Scrape down the sides of blender jar and break up the caked material with a spatula.



5. Add 150 ml petroleum ether and blend at high speed for 2 min.



6. Decant petroleum ether supernatant into a glassi fritted funnel with filter paper, into a 500 ml flask fitted with a suction apparatus.

7. Extract residues in blender with two additional 100 mL portions of petroleum ether, blending 2 min each time.
8. Decant petroleum ether supernatants through fritted funnel and combine with the first extract.
9. After last blending, transfer residue from blender jar to funnel, rinsing blender jar and materials in funnel with several small portions of petroleum ether.


10. Pour combined extracts and rinses into a column (25 mm x 50 mm long) of anhydrous  $\text{Na}_2\text{SO}_4$  and collect the petroleum ether extract.



11. Evaporate most of petroleum ether in a Kuderna-Danish concentrator.

12. Transfer the fat solution to tared pear-shaped flask using small amount of petroleum ether.

13. Evaporate petroleum ether at steam bath temperature using a rotary evaporator to obtain fat.



14. When petroleum ether is completely removed, weigh and record weight of fat extracted.


15. Take about 3 g fat for ether-acetonitrile partitioning.

16. Record weight of fat taken for clean-up using the formula:

$$\text{Wt. sample analyzed} = \frac{(\text{wt. fat for clean-up})}{(\text{wt. fat extracted})} \times (\text{wt. original sample})$$

## Clean-up

### a.) Petroleum ether-Acetonitrile Partitioning

1. Weigh about 3 g fat into a 125 mL separatory funnel.
  2. Add 15 mL petroleum ether and 30 mL acetonitrile saturated with petroleum ether.
  3. Shake the funnel vigorously for 1 min and allow the layers to separate.
- 
- The image shows a laboratory setting. A large, clear glass separatory funnel with a blue stopper and a yellow label '500 ml' is being tilted to pour a bright yellow liquid into a smaller, clear glass separatory funnel. The smaller funnel is positioned on a white lab bench. In the background, there are other lab equipment, including a beaker and some tubing.
4. Drain the acetonitrile into a 1 L separatory funnel containing 650 mL distilled H<sub>2</sub>O, 40 mL saturated NaCl solution, and 100 mL petroleum ether. Set aside.
  5. Extract the petroleum ether layer in 125 mL separatory funnel with 3 additional 30 mL portions of acetonitrile saturated with petroleum ether, shaking vigorously 1 min each time and combine all extracts in the 1L separatory funnel.
  6. Hold the 1 L separatory funnel in horizontal position and mix thoroughly for 30-45 sec. Allow the layers to separate and drain the aqueous layer into another 1 L separatory funnel.
  7. Add 100 mL petroleum ether to the second 1 L separatory funnel, shake vigorously for 15 sec and allow the layers to separate. Discard the aqueous layer, combine petroleum ether extract with petroleum ether in the first separatory funnel and wash with two 100 mL portions of water. Discard the washing.




8. Drain the petroleum ether layer through a column (25 mm x 50 mm long) of anhydrous  $\text{Na}_2\text{SO}_4$  into a 500 mL Kuderna-Danish concentrator.
9. Evaporate the extract to 5-10 mL in Kuderna-Danish concentrator and transfer the extract to a Florisil column.

#### **b.) Florisil Column Clean-up**

1. Prepare a column (10 mm i. d. x 120 mm long) and fill it with 4 g activated Florisil topped with anhydrous  $\text{Na}_2\text{SO}_4$  to about 2 cm. Pre-wet the column with 40-50 mL petroleum ether.







2. Transfer the petroleum ether solution of sample extract to the column. Rinse the container with two, 5 mL portions of petroleum ether and transfer rinsings to column.
3. Elute the column at about 5 mL/min with 200 mL of 6% ethyl ether/petroleum ether eluant.
4. Change the receiver and elute at about 5 mL/min with 200 mL of 15% ethyl ether/petroleum ether eluant.
5. Do step #4 again, this time using 50% ethyl ether/petroleum ether eluant.
6. Concentrate each eluate to a suitable definite volume (1-2 mL) in a Kuderna-Danish concentrator.
7. The concentrated eluate is now ready for injection into a gas chromatograph.

Notes: The first eluate (6%) contains organochlorine pesticides (aldrin, BHC, DDE, DDD, o,p'- and p,p'-DDT, heptachlor, heptachlor epoxide, lindane, methoxychlor, mirex, and ethylan), industrial chemicals (polychlorinated biphenyls (PCB) and organophosphorus pesticides (ethion and ronnel) and is usually suitable for gas chromatography directly.


The second eluate (15%) contains organochlorine pesticides (dieldrin and endrin) and organophosphorus pesticides (diazinon, methyl parathion and parathion).

The third eluate (50%) contains organophosphorus pesticide (malathion).

### c.) Magnesia Column Clean-up

(Applicable only to organochlorine pesticides in 15% eluate when additional clean-up is necessary).

1. Transfer about 10 g MgO celite mixture to a chromatographic tube without stopcock and using the vacuum to pack.



2. Pre-wash the column with about 40 mL petroleum ether, discard pre-wash, and place a Kuderna-Danish concentrator under the column.
3. Transfer the 15% Florisil eluate which has been concentrated to about 5 mL, to column, rinsing with small portions of petroleum ether.
4. Elute with 100 mL petroleum ether.
5. Concentrate the eluate to a suitable volume (1-2 mL). This is now ready for injection into a gas chromatograph.

## Detection and Quantitation by Gas Chromatography

Detection method is the process of identifying the pesticides and determining their concentrations. With the aid of an analytical instrument, such as a gas chromatograph (GC), we can obtain information about the original sample by running a standard solution containing the pesticide(s) at a known concentration and comparing it with the sample. In order for a sample to run through the GC, the pesticide must be extracted from the sample into a high purity solvent. This extracted portion is referred to as the sample extract. Upon entering the GC, the extract is vaporized in a thermally controlled injection port. The vapor containing the pesticide then passes into a temperature-and pressure-controlled oven within the GC where it is separated in the column. The vapor enters the column and is eventually transported, with the help of a carrier gas, to a thermally controlled detector which detects the pesticide as it flows through it. The time it takes for the pesticide to move through a column and reach the detector is its retention time and the same pesticide will always have the same retention time on the same column. All data transmitted from the detectors is collected and stored on a computer which compiles the data into a chromatogram as the GC is running.



## Gas Chromatographic Conditions

Gas chromatograph Model: Shimadzu GC-17A, equipped with  $^{63}\text{Ni}$ -electron capture detector, attached to a CBM-102 Chromatopak recorder system

Detector: Electron capture detector (ECD)

Column: SPB-608 (Supelco), Capillary (Fused silica), 30m x 0.25mm I. D., 0.25  $\mu\text{m}$  film

Column Oven Temperature: 150° C (4 min) to 290° C at 8° C/min, hold 10 min.

Detector temperature: 300° C

Injector temperature: 220° C

Carrier gas: Nitrogen ( $\text{N}_2$ ) at flow rate of 40 cm/sec

## Procedure

1. Study the operating manual of the gas chromatograph available in your laboratory. Turn on and set the specified chromatographic conditions.
2. Inject a suitable aliquot (3-8  $\mu\text{L}$ ) of concentrated eluate from Florisil or MgO-celite column containing an amount of compound within the linear range into the gas chromatograph.
3. Tentatively identify the residue peaks on basis of retention times.
4. Measure the area or height of residue peak(s) and determine the residue amount by comparison to peak area or height obtained from a known amount of appropriate reference material(s).
5. To ensure valid measurement of residue amount, the size of peaks from residue and reference standard should be within  $\pm 25\%$ .
6. Chromatograph reference material(s) immediately after samples.

## Determination of pesticide residues in non-fatty samples

### Principle

Non-fatty samples are blended with acetone and filtered. Extracts are transferred from aqueous filtrate to organic phase by shaking with petroleum ether and methylene chloride.

### Reagents

- a.) Solvents: acetone, methylene chloride ( $\text{CH}_2\text{Cl}_2$ ), petroleum ether- redistilled in glass
- b.) Sodium sulfate- Anhydrous, granular
- c.) Glass wool- Rinse with acetone and alcohol several times and dry.  
Washed glass wool will be somewhat brittle.
- d.) Reference standards- Prepare all standards in glass-distilled acetone

### Extraction

1. Weigh 100 g chopped sample into high-speed blender jar and add 200 mL acetone. Blend for 2 min at high speed.
2. Filter with suction through a Buchner funnel fitted with filter paper and collect the extract in 500 mL suction flask.
3. Place 80 mL sample extract in 1 L separatory funnel, and add 100 mL petroleum ether and 100 mL  $\text{CH}_2\text{Cl}_2$ . Shake vigorously for 1 min and let stand to allow layers to separate.
4. Transfer lower aqueous layer to a second 1 L separatory funnel. Set aside.
5. Dry the upper organic layer in first separatory funnel by passing through 1 ½ inch column of  $\text{Na}_2\text{SO}_4$  supported on washed glass wool in 4 inch funnel, and collect the eluate in 500 mL Kuderna-Danish concentrator fitted with volumetric flask or calibrated receiving tube.



6. To the aqueous phase in Step 4, add 7 g NaCl and shake vigorously until most NaCl is dissolved. Add 100 mL CH<sub>2</sub>Cl<sub>2</sub>, shake for 1 min, and dry lower organic phase through same Na<sub>2</sub>SO<sub>4</sub> column.
7. Extract the aqueous phase again with additional 100 mL CH<sub>2</sub>Cl<sub>2</sub> and dry as above. Rinse Na<sub>2</sub>SO<sub>4</sub> with 50 mL CH<sub>2</sub>Cl<sub>2</sub>.
8. Attach a Snyder column on Kuderna-Danish concentrator (boiling chips may be added) and start evaporation slowly by placing only receiver tube into steam.
9. After 100-150 mL has evaporated, concentrator may be exposed to more steam and concentrated further to about 2 mL.
10. Add 50 mL petroleum ether and repeat concentration step.
11. Add 20 mL acetone and reconcentrate to about 2 mL.

Note: Do not let solution go to dryness during any concentration step.

12. Adjust the volume of extract to suitable definite volume with acetone.
13. Inject a measured amount into the GC.

### Calculation of Equivalent Sample Weight

$$\frac{\text{mg sample equivalent}}{\text{uL final extract}} = 100 \times \frac{80}{200 + W - 10} \times \frac{1}{\text{mL final volume}}$$

where: 200 = mL acetone blended with 100 g sample

W = amount (mL) water present in original sample

10 = adjustment for water-acetone volume contraction

Note: Thus, when sample contains 85% water (85ml/100g) and final extract volume is 7 mL, each uL contains:

$$100 \times \frac{80}{200 + (85-10)} \times \frac{1}{7} = 4.15 \text{ mg sample equivalent/uL final extract}$$

## Determination of Polychlorinated Biphenyl Residues

### Extraction

1. Weigh 20 g of thoroughly ground and mixed sample.
2. Moisten 40 g granular  $\text{Na}_2\text{SO}_4$  with petroleum ether and add to the sample. Mix and let stand for 20 min and mix again.
3. Add 100 mL petroleum ether to the sample and blend for 1-2 min. Let stand for 10 min.
4. Plug funnel with glass wool, overlay with 20 g granular  $\text{Na}_2\text{SO}_4$ , and place funnel in 250 mL volumetric flask. Decant petroleum ether extract through  $\text{Na}_2\text{SO}_4$  into the volumetric flask.
5. Repeat the extraction of the sample using 100 mL petroleum ether.
6. Repeat extraction using 70 mL petroleum ether.
7. Dilute the filtrate by the addition of petroleum ether so that the total final volume is 250 mL.
8. Transfer 25 mL aliquot to a tared pear-shaped flask and evaporate the solvent in a rotary evaporator.
9. Weigh the flask and determine % fat.
10. For fish containing <10% fat, transfer a 25 mL aliquot to 125 mL K-D concentrator.  
For fish containing > 10% fat, take an aliquot to contain not > 200 mg fat.
11. Concentrate the sample to about 3 mL on a steam bath and transfer to a Florisil column.

### Clean-up (Florisil Column)

1. Place 4 g Florisil into (120 x 10 mm id) chromatographic tube. Add anhydrous  $\text{Na}_2\text{SO}_4$  to a height of 2 cm above the Florisil.
2. Open the stopcock fully, tap the tube to allow the adsorbent to settle, and mark the tube at 1 cm above the  $\text{Na}_2\text{SO}_4$  layer.

3. Wash the Florisil column with 20-25 mL petroleum ether.  
(Note: Solvent level must not go down below the 1-cm mark.)
4. Place 125 mL K-D flask under the column.
5. Using a disposable Pasteur pipet, transfer 3 mL sample to column. Temporarily close stopcock if necessary.
6. Add 35 mL petroleum ether-ether mixture (94 + 6) and elute.
7. Concentrate the eluate to about 2 mL.
8. Inject measured amount into GC.

## Determination of Carbamate Residues

### Principle

The carbamate group of pesticide residues includes carbonolates, carbaryl (Sevin), carbofuran, and propoxur (Baygon). The residue is extracted from the sample with acetonitrile. The extract is purified by partitioning with petroleum ether and coagulating in a  $\text{H}_3\text{PO}_4\text{-NH}_4\text{Cl}$  solution. Phenolic impurities are largely eliminated by partitioning  $\text{CH}_2\text{Cl}_2$  extract with KOH solution. Carbamate residues are treated with 1-fluoro-2,4-dinitrobenzene to form their corresponding derivatives.

### Reagents

- a.) Borax - 5% aqueous solution. Dissolve 5 g in 95 mL water.
- b.) Diatomaceous earth - Wash thoroughly with acetone and dry for 2 h at 110° C.
- c.) Coagulating solution:  
Stock solution- Dissolve 20 g  $\text{NH}_4\text{Cl}$  and 40 mL  $\text{H}_3\text{PO}_4$  in 360 mL distilled water.  
Working solution- Dilute 100 mL stock solution to 1 L for coagulation.
- d.) 1-fluoro-2,4 - dinitrobenzene solution - Redistill at 128° C and 1 mm pressure. Dissolve 1.5 mL in 25 mL acetone.

- e.) Pesticide standards - Best quality obtained from manufacturer, analytical grade when available
- f.) Potassium hydroxide solution - 0.5N aqueous solution
- g.) Sodium chloride solution - 30% aqueous solution
- h.) Solvents - acetone, methylene chloride, isooctane, acetonitrile and petroleum ether distilled in glass, acetophenone, and methanol (analytical grade)

### **Extraction and Clean-up**

1. Weigh 100 g sample in a beaker and add 200 mL  $\text{CH}_3\text{CN}$ .
2. Homogenize in blender operated at moderate speed for 2 min.
3. Filter with suction into 500 mL round-bottom flask through filter paper in 12 cm Buchner funnel.
4. Transfer 100 mL aliquot into 250 mL separatory funnel.
5. Add 25 mL NaCl solution and shake.
6. Drain and discard the aqueous phase and repeat the treatment with fresh NaCl solution.
7. Add 100 mL petroleum ether, and shake for 30 sec.
8. Drain  $\text{CH}_3\text{CN}$  into 1 L separatory funnel.
9. Strip petroleum ether by shaking for 20 sec with 50 and 10 mL portions of  $\text{CH}_3\text{CN}$ , draining each into the 1 L separatory funnel.
10. Add 300 mL distilled water, 25 mL NaCl solution, and 50 mL methanol.
11. Extract the mixture with 100 mL and two 25 mL portions of  $\text{CH}_2\text{Cl}_2$ , shaking each for 20 sec.
12. Drain the lower layer into 500 mL round-bottom flask.
13. Add 2 drops acetophenone, and evaporate in rotary evaporator under reduced pressure.

Note: During evaporation keep water bath within 40-50° C range and remove flask from water bath when the volume of the extract has been reduced to a few ml, so that final evaporation to dryness takes place at low temp.

14. Add 5 mL acetone, and swirl the flask to dissolve residue.
15. Add 50 mL coagulating solution, and swirl to mix.
16. Add 1-2 g diatomaceous earth, and swirl again to mix.
17. Pour the solution into a 150 mL suction filter of medium porosity and collect the filtrate in a 500 mL round-bottom flask.
18. Rinse the sides of flask with 5 mL acetone, swirl and repeat the coagulation.
19. When the filtration is completed, transfer the filtrate to 250 mL separatory funnel.
20. Extract the carbamate residue by shaking for 20 sec with three 25 mL portions of  $\text{CH}_2\text{Cl}_2$ .
21. Drain the  $\text{CH}_2\text{Cl}_2$  (lower layer) extract into another 250 mL separatory funnel. The solution may be kept overnight at this point.
22. Add 40 mL distilled water and 10 ml 0.5 N KOH, mix by gently and briefly swirling the flask, and shake for 20 sec.
23. Drain the  $\text{CH}_2\text{Cl}_2$  layer through granular anhydrous  $\text{Na}_2\text{SO}_4$  supported by glass wool in a filter funnel, and collect the filtrate in a 250 mL Erlenmeyer flask.
24. Add 100 mL distilled water, 2 mL 0.5 N KOH and 1 mL 1-fluoro-2,4-dinitrobenzene solution. Stopper, and mix for 20 min at high speed on a mechanical agitator.
25. Add 10 mL 5% borax, swirl to mix, and heat on a steam bath for 20 min.
26. Cool to room temperature by placing the flask in shallow water bath for 10 min.
27. Add 5 mL isooctane, stopper, shake for 3 min at high speed, and pour into 250 mL separatory funnel.
28. Drain the aqueous phase and rinse twice with distilled water.
29. Drain the isooctane solution into a funnel containing a 6 mm glass wool plug and into a test tube. Cover the test tube with a glass stopper. The solution may be kept overnight at this point.
30. Inject 10  $\mu\text{L}$  sample into a gas chromatograph. If it is necessary to dilute the sample, transfer 1 ml of the isooctane extract to another test tube, dilute to exact volume with isooctane, and shake to mix.



## Method Validation

The recovery of each pesticide residue was determined by spiking tiger shrimp or fish muscle tissue with 0, 5, 10, 20, 30, 40, 50, 100, 150 and 200 ng/g of each pesticide standard, and extracting the samples as described previously. Each determination was replicated 5 times. The recoveries ranged from 75 to 99%. Precision, expressed as percentage relative standard deviation, was below 2.5%.

### Limits of Detection

A method detection limit (MDL) was determined by running and analyzing a series of quality control samples (EPA STANDARDS, Spikes) to determine the lowest concentration of the pesticide that can be consistently and reliably measured using a given instrument. A spike is a standard containing the pesticide at a known concentration and goes through the entire extraction process like a real sample. Since spikes mimic real samples, we can measure the efficiency of the extraction method as well as the instrument performance and establish the (MDL) of the pesticide in various sample media. Laboratories should not report a pesticide at a concentration lower than the method detection limit.

An instrument detection limit (IDL) is determined by running a series of low level standards. The IDL gives us an idea of the lowest concentration of the pesticide residue that can be detected by the instrument. The IDL is consistently lower than the MDL because there is no extraction involved. We only use IDLs to monitor the instrument, but the IDLs have little relevance to the concentration of the actual sample extract.

The MDLs of 29 pesticide residues established using the methods described previously are summarized in the table next page:

## Method Detection Limit (MDL) of 29 Pesticide Residues

Pesticides	MDL µg/g	Pesticides	MDL µg/g
Aldrin	0.075	Endrin ketone	0.024
a- BHC	0.025	Heptachlor	0.010
b- BHC	0.010	Heptachlor epoxide	0.015
g- BHC	0.010	Methoxychlor	0.050
d- BHC	0.015	Chloroneb	0.50
a-Chlordane	0.0015	Chlorothanonil	0.025
g- Chlordane	0.0015	Etridiazole	0.025
p,p' DDD	0.0025	Propachlor	0.050
p,p' DDE	0.010	Trifluralin	0.025
p,p' DDT	0.060	Simazine	0.50
Dieldrin	0.020	Aldicarb	1.0
Endosulfan I	0.015	Carbaryl (Sevin)	2.0
Endosulfan II	0.015	Carbofuran	1.5
Endosulfan sulfate	0.015	Propour (Baygon)	1.0
Endrin	0.025		

## References

- (AOAC). Official Methods of Analysis. 1996. Sawyer, L. D., MacMahon, B. M., and Newsome, W. H. (Chapter eds.). Pesticides and Industrial Chemical Residues. 10:1-78. AOAC International, 16<sup>th</sup> edition. Arlington, VA, U.S. A.
- Bertuzzi, P. F., Kamps, L. and Miles, C. I. 1967. Extraction of chlorinated pesticide residues from non-fatty samples of low moisture content. J. Assoc. Offic. Agr. Chemists. 50:623-627.
- Cahill, W. P., Estes, B. J., and Ware, G. W. 1970. A rapid on-column extraction-clean-up method for animal fat. Bull. Environ. Contam. Toxicol. 5: 70-71.
- Carr, R. L. 1971. Collaborative study of a method for multiple chlorinated pesticide residues in fish. J. Assoc. Offic. Agr. Chemists. 54:525-527.
- Jain, C. K. and Ali, I. 1997. Determination of pesticides in water, sediments, and soils by gas chromatography. Int. J. Environ. Anal. Chem. 68(1): 83-101.
- Porter, M. L., Young, S. V. and Burke, J. A. 1970. A method for the analysis of fish, animal and poultry tissue for chlorinated pesticide residues. J. Assoc. Offic. Agr. Chemists. 50:644-645.